

INTOXILYZER 8000

Location: JAIL

Serial Number: 80-003355

Core Version: 8105.48

08/22/2018

22:03:58

Standard Lot#: 13917100A1

Last Changed By: D. SPIVEY #T17

QAS: J. SONDROL #187

LAKE HAVASU CITY PD

Operator: J. SONDROL #187

LAKE HAVASU CITY PD

Test	g/210L	Time
Air Blank	0.000	22:04:30
0.100 Cal Check	0.100	22:04:51
Air Blank	0.000	22:05:21

Preventative Maintenance
Performed -- Timer Reset

Exhibit G-2. Standard Calibration Check Procedure, Intoxilyzer Model 8000

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J SONDROL 187 AGENCY LHC PD
 DATE 8/22/18 TIME 2203
 INTOXILYZER SERIAL # 80-003355 LOCATION JAIL

() 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
 OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34° C ± 0.2° C

- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
 () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
 () 4. Air blank completed.
 () 5. Calibration check completed. Test results 0. 100 AC.
 () 6. Air blank completed.
 () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J Sondrol 187

DPS Form Exh G-2 (Rev 05-01)

Historical Note

New Exhibit G-2 made by final rulemaking at 12 A.A.R. 1916, effective 9:00 a.m., May 18, 2006 (Supp. 06-2).