



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MALOY #281 AGENCY LKCPD  
DATE 11-24-2021 TIME 1855  
INTOXILYZER SERIAL # 80-003355 LOCATION JAIL

- (  ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (  ) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (  ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (  ) 4. Air blank completed.
- (  ) 5. Calibration check completed. Test results 0. 097 AC.
- (  ) 6. Air blank completed.
- (  ) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Michael Maloy #281

INTOXILYZER 8000  
Location: JAIL  
Serial Number: 80-003355  
Core Version: 8105.48  
11/24/2021 18:55:15

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281

-----  
GAS: M. MALOY JR #281  
LAKE HAVASU CITY PD

-----  
Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD

-----

Test	g/210L	Time
Air Blank	0.000	18:55:47
0.100 Cal Check	0.097	18:56:08
Air Blank	0.000	18:56:37

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. MALOY #281 AGENCY LKCPA  
DATE 11-24-2021 TIME 1859  
INTOXILYZER SERIAL # 80-003355 LOCATION JAIL

(  ) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (  ) 1. Clock time check.
- (  ) 2. Date check.

OPERATIONAL TESTS

- (  ) 1. Alcohol-free subject test result 0. 000 AC.
- (  ) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- (  ) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- (  ) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

RETURNED INSTRUMENT BACK INTO SERVICE FROM SCHEDULED P.M.

SIGNATURE Michael Maloy Jr #281

INTOXILYZER 8000  
 Location: JAIL  
 Serial Number: 80-003355  
 Core Version: 8105.48  
 11/24/2021 18:59:02

Standard Lot#: 35019100A1  
 Last Changed By: M. MALOY JR #281

OAS: M. MALOY JR #281  
 LAKE HAVASU CITY PD

Operator: M. MALOY JR #281  
 LAKE HAVASU CITY PD

Subject: SOAP,TEST  
 DOB: 01/01/2000  
 Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:59:50
Diagnostic Test	Pass	19:00:28
Air Blank	0.000	19:00:55
0.100 Cal Check	0.096	19:01:16
Air Blank	0.000	19:01:46
Subject Test	0.000	19:02:07
Air Blank	ABT*	19:02:13
Air Blank	1PS**	19:02:21

\*Sequence Aborted  
 \*\*Improper Sample

Not a Successfully  
 Completed  
 Test Sequence

INTOXILYZER 8000  
 Location: JAIL  
 Serial Number: 80-003355  
 Core Version: 8105.48  
 11/24/2021 19:05:20

Standard Lot#: 35019100A1  
 Last Changed By: M. MALOY JR #281

OAS: M. MALOY JR #281  
 LAKE HAVASU CITY PD

Operator: M. MALOY JR #281  
 LAKE HAVASU CITY PD

Subject: SOAP,TEST  
 DOB: 01/01/2000  
 Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:06:12
Diagnostic Test	Pass	19:06:50
Air Blank	0.000	19:07:17
0.100 Cal Check	0.097	19:07:38
Air Blank	0.000	19:08:08
Subject Test	0.000*	19:11:15
Air Blank	0.000	19:11:43
Air Blank	0.000	19:12:10
Subject Test	XXX**	19:12:23
Air Blank	RFI***	19:12:35

\*Deficient Sample  
 \*\*Invalid Sample  
 Begin new deprivation period  
 \*\*\*RFI Detect

Not a Successfully  
 Completed  
 Test Sequence

281

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MALOY #281 AGENCY LNCPD  
DATE 11-24-2021 TIME 1912  
INTOXILYZER SERIAL # 80-003355 LOCATION JAIL

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Michael Maloy #281

INTOXILYZER 8000  
Location: JAIL  
Serial Number: 80-003355  
Core Version: 8105.48  
11/24/2021 19:12:56

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281  
-----

OAS: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Test	g/210L	Time
Air Blank	0.000	19:13:28
0.100 Cal Check	0.096	19:13:49
Air Blank	0.000	19:14:19

Preventative Maintenance  
Performed -- Timer Reset