

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. S. [unclear] 157 AGENCY LW 83

DATE 7/20/17 TIME 5:00

INTOXILYZER SERIAL # 80-059336 LOCATION MAINE JAIL

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 0.08 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed check list.

SIGNATURE [Signature]

4/11

INTOXILYMER 8000

Location: MAIN JAIL

Serial Number: 80-003156

Core Version: 9105.48

10/13/2017

02:42:18

Standard ID#: 10318 0001

Last Checked By: J. SONDRON #187

QAS: J. SONDRON #187

LAKE HAVASU CITY PD

Operator: J. SONDRON #187

LAKE HAVASU CITY PD

Line	Time	Time
01	11:17:50	
02	11:18:09	
03	11:18:48	

Printed at the Maintenance

10/13/2017 02:42:18

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. Salazar 137 AGENCY LNC PD
DATE 7/24/15 TIME 2044
INTOXILYZER SERIAL # 80003354 LOCATION PT. RAIN STATION

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0.000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0.111 AC. Result 0.111 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE J. Salazar

INTOXILYZER 8000
Location: MAIN JAIL
Serial Number: 80-003356
Core Version: 8105.48
07/20/2017

22:44:31

Standard Lot#: 22315100A1
Last Changed By: J. SONDROL #187

QAS: J. SONDROL #187
LAKE HAVASU CITY PD

Operator: J. SONDROL #187
LAKE HAVASU CITY PD

Subject: SQAP, TEST

DOB: 01/01/1950

Sex: M

Weight: 123

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:46:27
Diagnosed Test	Pass	22:47:03
Air Blank	0.000	22:47:29
0.100 Cal Check	0.097	22:47:56
Air Blank	0.000	22:48:20
Sampled Test	0.000	22:48:49
Air Blank	APF*	22:49:03
Air Blank	APF*	22:49:11

Acquisition Aborted

Sampled Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: MAIN JAIL

Serial Number: 80-003356

Core Version: 8105.48

07/20/2017

22:49:54

Standard Lot#: 22315100A1

Last Changed By: J. SONDRON #187

QAS: J. SONDRON #187

LAKE HAVASU CITY PD

Operator: J. SONDRON #187

LAKE HAVASU CITY PD

Subject: SQAP, TEST

DOB: 01/01/1950

Sex: M

Weight: 123

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:50:38
Diagnostic Test	Pass	22:51:13
Air Blank	0.000	22:51:40
0.100 Cal Check	0.097	22:52:01
Air Blank	0.000	22:52:31
Subject Test	0.175 ^U	22:55:38
Air Blank	0.000	22:56:12
Air Blank	0.000	22:56:39
Subject Test	xxxx ^U	22:57:03
Air Blank	0.000	22:57:57

**Inefficient Sample

**Inverted Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. J. Williams AGENCY Phoenix

DATE 12/11 TIME 3:00

INTOXILYZER SERIAL # 30-00000 LOCATION 118

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. _____ AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: MAIN JAIL

Serial Number: 80-003356

Core Version: 8105.48

07/20/2017

23:00:48

Standard Lot#: 22315'00A'

Last Changed By: J. SONDRON #187

QAS: J. SONDRON #187

LAKE HAVASU CITY PD

Operator: J. SONDRON #187

LAKE HAVASU CITY PD

Test	g 210L	Time
Alt Blank	0.000	23:01:17
0.100 211 Check	0.100	23:01:38
Alt Blank	0.000	23:02:07

Preventative Maintenance

Performed: Check 1-1-17