



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MALDY #281 AGENCY LUCPD

DATE 03-27-2021 TIME 0943

INTOXILYZER SERIAL # 80-001873 LOCATION STATION

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 097 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Michael Maloy Jr #281

281

INTOXILYZER 8000  
Location: LACPD  
Serial Number: 80-001873  
Core Version: 8105.48  
03/27/2021 09:43:12

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281  
-----

QAS: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Test	g/210L	Time
Air Blank	0.000	09:43:43
0.100 Cal Check	0.097	09:44:04
Air Blank	0.000	09:44:33

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. MALDY #281 AGENCY LHCPO

DATE 03-27-2021 TIME 0949

INTOXILYZER SERIAL # 80-001873 LOCATION STATION

(  ) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(  ) 1. Clock time check.

(  ) 2. Date check.

OPERATIONAL TESTS

(  ) 1. Alcohol-free subject test result 0. 000 AC.

(  ) 2. Error recognition logic system functioning.

Not a Successfully Completed Test Sequence printed

(  ) 3. Proper sample recognition system.

Not a Successfully Completed Test Sequence printed

Deficient sample printed.

(  ) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

RETURNED TO AGENCY FROM AZ O.P.S. IN  
PHOENIX FOR REPAIR AND PREVENTATIVE  
MAINTENANCE.

SIGNATURE Michael Maloy Jr

INTOXILYZER 8000  
Location: LHCPD  
Serial Number: 80-001873  
Core Version: 8105.48  
03/27/2021 09:49:15

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281  
-----

-----  
OAS: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

-----  
Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

-----  
Subject: SQAP,TEST  
DOB: 01/01/2000  
Sex: M Weight: 200  
-----

-----  
15 Minute Deprivation Period? Yes  
-----

Test	g/210L	Time
Air Blank	0.000	09:50:04
Diagnostic Test	Pass	09:50:42
Air Blank	0.000	09:51:09
0.100 Cal Check	0.097	09:51:30
Air Blank	0.000	09:52:00
Subject Test	0.000	09:52:16
Air Blank	ABT*	09:52:23
Air Blank	IPS**	09:52:30

\*Sequence Aborted  
\*\*Improper Sample

Not a Successfully  
Completed  
Test Sequence

INTOXILYZER 8000  
Location: LHCPD  
Serial Number: 80-001873  
Core Version: 8105.48  
03/27/2021 09:59:14

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281  
-----

-----  
OAS: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

-----  
Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

-----  
Subject: SQAP,TEST  
DOB: 01/01/2000  
Sex: M Weight: 200  
-----

-----  
15 Minute Deprivation Period? Yes  
-----

Test	g/210L	Time
Air Blank	0.000	10:00:04
Diagnostic Test	Pass	10:00:42
Air Blank	0.000	10:01:09
0.100 Cal Check	0.097	10:01:30
Air Blank	0.000	10:02:00
Subject Test	0.000*	10:05:07
Air Blank	0.000	10:05:34
Air Blank	0.000	10:06:02
Subject Test	XXX**	10:06:16
Air Blank	RFI***	10:06:30

\*Deficient Sample  
\*\*Invalid Sample  
Begin new deprivation period  
\*\*\*RFI Detect

Not a Successfully  
Completed  
Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MALOY #281 AGENCY LHCPD

DATE 03-27-2021 TIME 1008

INTOXILYZER SERIAL # BD-001873 LOCATION STATION

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 097 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Michael Maloy #281

INTOXILYZER 8000  
Location: LHCPD  
Serial Number: 80-001873  
Core Version: 8105.48  
03/27/2021 10:08:03

-----  
Standard Lot#: 35019100A1  
Last Changed By: M. MALOY JR #281  
-----

QAS: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Operator: M. MALOY JR #281  
LAKE HAVASU CITY PD  
-----

Test	g/210L	Time
Air Blank	0.000	10:08:37
0.100 Cal Check	0.097	10:08:58
Air Blank	0.000	10:09:27

Preventative Maintenance  
Performed -- Timer Reset