

To Whom It May Concern:

I hereby attest that Larry Allen, is a Certified Intoxilyzer Quality Assurance Specialist for the Lake Havasu City Police Department.

I certify that the attached documentation of Standard Quality Assurance Procedures and/or Standard Calibration Check Procedure are full, true and correct copy of the original.

Signature

COUNTY OF MOHAVE) ss)	
Subscribed and sworn to before 7020 by 440	ore me this day of	MARCH.
NOTARY PUBLIC		
03-11-2020 DATE		

Commission expires:

STATE OF ARIZONA



Area Code 928

Administration	Fax #'s Administration	Emergency
Patrol 855-0515	Investigations/Patrol 680-5432	

E-mail: police@lhcaz.gov

EXHIBIT G-2 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST L. ALLEN #251 AGENCY LHC PA
DATE 03-11-2020 TIME 0033
INTOXILYZER SERIAL # 80-003355 LOCATION LAC PD JA; L
() 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution OR
Pour a standard alcohol concentration solutionAC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
 Intoxilyzer 8000 display reads "PUSH BUTTON TO START" Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
 4. Air blank completed. 5. Calibration check completed. Test results 0. 099 AC.
6. Air blank completed. 7. Remove printed record. Attach the record to the completed checklist.
SIGNATURE # Z51
DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000 Location: JAIL

Serial Number: 80-003355 Core Version: 8105.48

03/11/2020 00:33:59

Standard Lot#: 13518100A1

Last Changed By: J. SONDROL #187

OAC. T. NIEN #251

QAS: L. ALLEN #251 LAKE HAVASU CITY PD

Operator: L. ALLEN #251 LAKE HAVASU CITY PD

.

Test	g/210L	Time
Air Blank	0.000	00:34:37
0.100 Cal Check	0.099	00:34:58
Air Blank	0.000	00:35:27

Preventative Maintenance
Performed -- Timer Reset

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST L ALLEN # ZSI AGENCY LHC PA
DATE 03-11-2020 TIME 0042
INTOXILYZER SERIAL # 80-003355 LOCATION LAC PD JAIL
(/) 1. Display Reads "PUSH BUTTON TO START"
DIAGNOSTIC TESTS () 1. Clock time check. () 2. Date check.
OPERATIONAL TESTS (
Instrument is operating properly and accurately. YES NO
COMMENTS
SIGNATURE #Z51
DPS Form Exh G-4 (Rev 05-01)

ケンン

INTOXILYZER 8000 Location: JAIL

Serial Number: 80-003355 Core Version: 8105.48

03/11/2020

00:42:08

Standard Lot#: 13518100A1 Last Changed By: J. SONDROL #187 QAS: L. ALLEN #251 LAKE HAVASU CITY PD Operator: L. ALLEN

LAKE HAVASU CITY PD

#251

Subject: SQAP, TEST DOB: 01/01/1995

Sex: M

Weight: 170 ·

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:44:02
Diagnostic Test	Pass	00:44:38
Air Blank	0.000	00:45:05
0.100 Cal Check	0.099	00:45:26
Air Blank	0.000	00:45:55
Subject Test	0.000	00:46:15
Air Blank	ABT*	00:46:26
Air Blank	IPS**	00:46:35

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

^{**}Improper Sample

INTOXILYZER 8000 Location: JAIL

Serial Number: 80-003355 Core Version: 8105.48

03/11/2020

00:48:01

Standard Lot#: 13518100A1

Last Changed By: J. SONDROL

#187

QAS: L. ALLEN #251

LAKE HAVASU CITY PD

.....

Operator: L. ALLEN

#251

LAKE HAVASU CITY PD
Subject: SQAP, TEST

DOB: 01/01/1995

Sex: M

Weight: 170

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:49:28
Diagnostic Test	Pass	00:50:03
Air Blank	0.000	00:50:30
0.100 Cal Check	0.099	. 00:50:51
Air Blank	0.000	00:51:21
Subject Test	0.000*	00:54:28
Air Blank	0.000	00:54:56
Air Blank	0.000	00:55:23
Subject Test	XXX**	00:55:45
Air Blank	RFI***	00:56:10

^{*}Deficient Sample

Begin new deprivation period

^{**}Invalid Sample

^{***}RFI Detect

EXHIBIT G-2 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST L ALLEW # 251 AGENCY LHCPD
DATE 03-11-2020 TIME 0057
INTOXILYZER SERIAL # 80-003355 LOCATION LHC PD JAIL
1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solutionAC
Pour a standard alcohol concentration solutionAC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
 Intoxilyzer 8000 display reads "PUSH BUTTON TO START" Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed. () 5. Calibration check completed. Test results 0. 098 AC. () 6. Air blank completed. () 7. Remove printed record. Attach the record to the completed checklist.
6. Air blank completed. 7. Remove printed record. Attach the record to the completed checklist.
SIGNATURE #751
DPS Form Exh G-2 (Rev 05-01)

#151

INTOXILYZER 8000 Location: JAIL

Serial Number: 80-003355 Core Version: 8105.48

03/11/2020 00:57:35

Standard Lot#: 13518100A1

Last Changed By: J. SONDROL #187

QAS: L. ALLEN #251 LAKE HAVASU CITY PD

Operator: L. ALLEN #251 LAKE HAVASU CITY PD

Test	g/210L	Time
Air Blank	0.000	00:58:10
0.100 Cal Check	0.098	00:58:31
Air Blank	0.000	00:59:01

Preventative Maintenance Performed -- Timer Reset