

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J Sondrol 187 AGENCY LNC PD
DATE 9/22/16 TIME 0025
INTOXILYZER SERIAL # 80-003355 LOCATION EO BOOKING

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

8187

SIGNATURE J Sondrol 187

DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000
Instrument Initialization
00:18:35 09/22/2016

INTOXILYZER 8000
Location: EAST BANK BOOKING
Serial Number: 80-003355
Concentration: 0.099 AC
09/22/2016 00:25:27

Standard Lot#: 22315100A1
Last Changed By: J. SONDROL #187

DPS: J. SONDROL #187
LAKE HAVASU CITY PD

Operator: J. SONDROL #187
LAKE HAVASU CITY PD

Test	g/210L	Time
Air Blank	0.000	00:25:56
0.100 Cal Check	0.099	00:26:17
Air Blank	0.000	00:26:46

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. SONOROL #187 AGENCY LNC PD
 DATE 9/22/16 TIME 0027
 INTOXILYZER SERIAL # 80-003355 LOCATION EE BOOKING

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
 () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
 () 2. Error recognition logic system functioning.
 Not a Successfully Completed Test Sequence printed
 () 3. Proper sample recognition system.
 Not a Successfully Completed Test Sequence printed
 Deficient sample printed.
 () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

g/s157

Instrument is operating properly and accurately. YES NO

COMMENTS

CRIME BACK FROM 90 DAY MAINTENANCE. COMPLETE
90-DAY

SIGNATURE *J. Sonorol*

INTOXILYZER 8000
 Location: EAST BANK BOOKING
 Serial Number: 80-003355
 Core Version: 8105.48
 09/22/2016 00:27:04

Standard Lot#: 22315100A1
 Last Changed By: J. SONOROL #187

DPS: J. SONOROL #187
 LAKE HAVASU CITY PD

Operator: J. SONOROL #187
 LAKE HAVASU CITY PD

Subject: SCAP, TEST
 DOB: 01/01/1950
 Sex: M Weight: 123

15 Minute Deprivation Period? Yes

DPS Form Exh G-4 (Rev 05-01)

Test	g/210L	Time
Air Blank	0.000	00:26:51
Diagnostic Test	Pass	00:26:53
Air Blank	0.000	00:26:56
1.00 Cal Check	0.999	00:26:56
Air Blank	0.000	00:27:26
Subject Test	0.000	00:27:26
Air Blank	0.000	00:27:46
Air Blank	0.000	00:27:46
Subject Test	0.000	00:27:46
Air Blank	0.000	00:27:46

Test	g/210L	Time
Air Blank	0.000	00:27:50
Diagnostic Test	Pass	00:28:27
Air Blank	0.000	00:28:55
1.00 Cal Check	0.999	00:29:15
Air Blank	0.000	00:29:46
Subject Test	0.000	00:30:04
Air Blank	0.000	00:31:22
Air Blank	0.000	00:31:30

*Sequence aborted
 **Inhibitor Sample

Deficient Sample
 Inhibitor Sample
 Best Test Method Used
 Inhibitor Sample

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J SONDROL 187 AGENCY LHC
DATE 9/22/16 TIME 0842
INTOXILYZER SERIAL # 80-003355 LOCATION EB BOOKING

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J Sondrol 187

DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000
Location: EAST BANK BOOKING
Serial Number: 80-003355
Core Version: 8105.48
09/22/2016 00:42:31

Standard Lot#: 22315100A1
Last Changed By: J. SONDROL #187

QAS: J. SONDROL #187
LAKE HAVASU CITY PD

Operator: J. SONDROL #187
LAKE HAVASU CITY PD

Test	g/210L	Time
Air Blank	0.000	00:43:00
0.100 Cal Check	0.097	00:43:21
Air Blank	0.000	00:43:50

Preventative Maintenance
Performed -- Timer Reset

